

SONOderma™

Aquabrasion™
Ultrasonic Skin Renewal



Revolutionary advance in the care of aging skin

brera
medical inc USA

AquabrasionTM using SONOdermaTM



SONOderma is a revolutionary advance in the care of aging skin. Using patented ultrasound technology and clean water to create a localized cavitation effect to gently break up the skin's corneous layer and unblock pores, this non-invasive system delicately exfoliates and rinses away dead skin cells, excess sebum, and bacteria, revealing younger, fresher skin that feels softer and cleaner.



SONOderma's particle-free Aquabrasion process is a unique, in-office dermabrasion procedure that can be performed by any trained, skin-care professional. Typically, treatments are performed up to six times every two to three weeks, and then every one to two months as part of a personalized anti-aging maintenance program.


Revolutionary Technology

- Advanced low-frequency ultrasonic transducer technology
- Multiple modes of action offer a range of aesthetic treatments
- Compact system, expandable with other technologies for greater versatility
- Can be used in conjunction with other facial products and treatments
- Safely and effectively exfoliate all skin types (including sensitive skin) without the use of abrasive particles
- Pleasant, painless treatments
- Immediate, in-demand results


Modes of Actions

- Gentle, mechanical exfoliation smooths and renews dull, dry skin
- Ultrasonic massage revitalizes skin with, or without the simultaneous use of topical skin-care products

Abrasion of corneous layer



Increasing cutaneous permeability



brera

medical inc USA sales@breramedicalusa.com

Aquabrasion™: Advancing Microdermabrasion

Microdermabrasion has been an important part of any practice dedicated to facial aesthetics. But now, SONOderma™ takes exfoliation technology to a new level using ultrasonic transducers instead of crystals and vacuum.

Some of the reasons why medical professionals prefer SONOderma Aquabrasion over traditional microdermabrasion include:

- No messy crystals
- No consumables
- No discomfort
- No deliberate skin irritation
- Can be used on sensitive skin, including rosacea and painful acne
- Can be used to gently massage the skin and in conjunction with a myriad of skin care products and treatments
- The Imperium® system is expandable to radio frequency (RF) and other ultrasound modalities



SONOderma is part of the **Imperium System** from **Brera® Medical**


The Imperium system is a compact, versatile and expandable aesthetic system used for many applications related to the rejuvenation aging skin. In addition to SONOderma, Imperium offers several different RF handpieces, as well as additional ultrasound handpieces that offer a wider range of treatment options, all related to the aesthetic care of the face and body. Please contact Brera Medical today for specific information regarding these additional modalities.




Imperium® Non-Invasive Face and Body System

Featuring
SONOderma™ Aquabrasion™
Ultrasonic Skin Renewal

Before



After



Specifications

Vibration

Frequency: 25,000 Hz ($\pm 10\%$)

Applicator: Medical grade stainless steel,
1 $\frac{3}{8}$ " x 1 $\frac{3}{4}$ " x $\frac{1}{2}$ " mm

Dimensions: 17.7" (W) x 13.4" (D) x 5.9" (H)

Weight: 19.8 lbs

Electrical: 110–230 VAC ($\pm 10\%$), 50/60 Hz

Imperium **brera**
medical inc USA

Brera Medical Inc. • 7 Upson Road, Wellesley, MA 02482 • 1 (888) 580-1233 • www.brerammedicalusa.com
© Brera Medical Inc, 2010. All rights reserved. Brera and Imperium are registered trademarks of Brera Medical Technologies.
SONOderma and Aquabrasion are trademarks of Brera Medical Inc.